Tiefziehtechnik
Thermoforming Technique
Technique de Thermoformage
Tecnica de termomoldeo
Tecnica di termoformatura

Material • Materials • Matériaux • Materiales • Materiale

	(D)	(GB)	(F)	(E)	(I)
BIOCRYL® C	4	16	28	40	52
BIOPLAST®	4	16	28	40	52
BIOPLAST® bleach	5	17	29	41	53
BIOPLAST® Color/Multicolor	5	17	29	41	53
BIOPLAST® XTREME	6	18	30	42	54
COPYPLAST®	6	18	30	42	54
COPYPLAST® C	7	19	31	43	55
DURAN®	7	19	31	43	55
DURAN®+	8	20	32	44	56
DURAN®+ A2 pd	8	20	32	44	56
DURAN®+ white pd	9	21	33	45	57
DURASOFT® pd	9	21	33	45	57
DURASOFT® blue pd	10	22	34	46	58
DURASOFT® seal	10	22	34	46	58
HARDCAST®	11	23	35	47	59
IMPRELON® natur	11	23	35	47	59
IMPRELON® klar	12	24	36	48	60
IMPRELON® S pd / IMPRELON® S+ pd	12	24	36	48	60
ISOFOLAN®	13	25	37	49	61
KASCHIERFOLIE	13	25	37	49	61
PLATZHALTERFOLIE	14	26	38	50	62
UNIQ pd / UNIQ+ pd	14	26	38	50	62

BIOCRYL® C

Indikation

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CE 0044

BIOCRYL® C**Indication**

Immediate dentures / temporary dentures and retention plate / expansion plate

Tips and tricks:

Use of ISOFOLAN® required, wire elements and screws can be polymerised using BIOCRYL®-RESIN, take care while polishing (brittle)

Available in:

2.0 / 3.0 mm x 125 mm Ø
clear

2.0 mm x 125 mm Ø
red, blue, yellow, pink-transparent

Technical data

Material: PMMA
E-modulus: 3,300 MPa
Hardness: Rockwell 195
Tensile strength: 70 MPa
Impact strength: 20 KJ/m²

C € 0044

BIOPLAST®**Indication**

Duplicating form for situation models and OSAMU-Retainer® and positioner and sports mouthguard

Tips and tricks:

Can be melted together using IR heater or heat gun. Thermal polishing possible. Cleaning and perforation required. Use of BIOPLAST® insulating agent recommended.

Available in:

1.0 / 1.5 / 2.0 / 3.0 / 4.0 / 5.0 mm
x 125 mm Ø

Technical data

Material: Ethylene-vinyl acetate
E-modulus: 15 MPa
Hardness: approx. 78–82 Sh. A
Tensile strength: 18 MPa
Impact strength: no fracture

C € 0044

BIOPLAST® bleach**Indication**

Bleaching splint

Tips and tricks:

Use of BIOPLAST® insulating agent recommended.

Available in:

1.0 mm x 125 mm Ø

Technical data

Material: Ethylene-vinyl acetate
E-modulus: 20 MPa
Hardness: approx. 92 Sh. A
Tensile strength: 20 MPa
Impact strength: no fracture

C €

BIOPLAST® Color/Multicolor**Indication**

Sports mouthguard

Tips and tricks:

Can be melted together using IR heater or heat gun. Thermal polishing possible. Cleaning and perforation required. Use of BIOPLAST® insulating agent recommended.

Available in:

3.0 mm x 125 mm Ø
Various colour combinations

Technical data

Material: Ethylene-vinyl acetate
E-modulus: 20 MPa
Hardness: approx. 78–82 Sh. A
Tensile strength: 18 MPa
Impact strength: no fracture

C €

BIOPLAST® XTREME

Indication

Sports mouthguard

Tips and tricks:

Additional reinforcement can be achieved by laminating DURAN® and BIOPLAST®. Cleaning and perforation required, use of BIOPLAST® Insulating agent recommended.

Available in:

4.0 / 5.0 mm x 125 mm Ø
red, blue, black and clear

Technical data

Material:	Ethylene-vinyl acetate
E-modulus:	15/20 MPa
Hardness:	80/90 Sh. A
Tensile strength:	18/20 MPa
Impact strength:	no fracture

COPYPLAST®

Indication

Etching tray and bleaching splint and bracket-transfer mask / transfer splint and duplicating form for situation models

Tips and tricks:

Insulation against all kinds of resins, also suitable for the production of PMMA duplicates.

Available in:

0.5 / 0.75 / 1.0 / 1.5 / 2.0 mm
x 125 mm Ø

Technical data

Material:	Polyethylene
E-modulus:	175 MPa
Hardness:	42 Sh. D
Tensile strength:	9 MPa
Impact strength:	101 KJ/m²

COPYPLAST® C

Indication

Etching tray and bleaching splint and retainer / retention splint / long-term retainer

Tips and tricks:

Can be trimmed with scissors or scalpel, edges can be smoothed by means of DIMO® wheel or silicone polisher.

Available in:

1.0 mm x 125 mm Ø

Technical data

Material:	Polypropylen
E-modulus:	450 MPa
Hardness:	68 Sh. D
Tensile strength:	19 MPa
Impact strength:	650 KJ/m²

DURAN®

Indication

Aligner and occlusal splint and fluorisation splint and interim splint and orthodontic correction splint and bruxism splint / miniplast splint and retainer / retention splint / long-term retainer and surgical X-ray splint and sports mouthguard and Mandibular Protrusion Splint / MAD / anti-snoring device

Tips and tricks:

Occlusal adjustments can be done using DURASPLINT® without bonding agent, use of ISOFOLAN® required.

Available in:

0.4 / 0.5 / 0.625 / 0.75 / 1.0 / 1.5 / 2.0 / 3.0 mm x 125 mm Ø

Technical data

Material:	PETG
E-modulus:	2,200 MPa
Hardness:	78 Sh. D
Tensile strength:	53 MPa
Impact strength:	10 KJ/m²

DURAN®+**Indication**

Aligner and occlusal splint and fluorisation splint and orthodontic correction splint and bruxism splint / miniplast splint and retainer / retention splint / long-term retainer and surgical X-ray splint and sports mouthguard and Mandibular Protrusion Splint / MAD / anti-snoring device

Tips and tricks:

Occlusal adjustments can be done using DURASPLINT® without bonding agent. Comes with integrated insulation foil that is removed after thermoforming, use of ISOFOLAN® not required.

Available in:

0.4 / 0.5 / 0.625 / 0.75 / 1.0 / 1.5 / 2.0 mm
x 125 mm Ø ; with integrated protection foil

Technical data

Material:	PETG
E-modulus:	2,200 MPa
Hardness:	78 Sh. D
Tensile strength:	53 MPa
Impact strength:	10 KJ/m ²

C € 0044

DURAN®+ A2 pd**Indication**

Interim splint and aesthetic splint

Tips and tricks:

Occlusal adjustments can be done using DURASPLINT® without bonding agent. Comes with integrated insulation foil that is removed after thermoforming, use of ISOFOLAN® not required.

Available in:

0.5 / 0.625 / 0.75 mm
x 125 mm Ø ; with integrated protection foil

Technical data

Material:	PETG
E-modulus:	2,200 MPa
Hardness:	78 Sh. D
Tensile strength:	53 MPa
Impact strength:	10 KJ/m ²

C € 0044

DURAN®+ white pd**Indication**

Aesthetic splint

Tips and tricks:

Occlusal adjustments can be done using DURASPLINT® without bonding agent, comes with integrated insulation foil that is removed after thermoforming, use of ISOFOLAN® not required.

Available in:

0.5 / 0.625 / 0.75 / 1.0 mm
x 125 mm Ø ; with integrated protection foil

Technical data

Material:	PETG
E-modulus:	2,200 MPa
Hardness:	78 Sh. D
Tensile strength:	53 MPa
Impact strength:	10 KJ/m ²

C €

DURASOFT® pd**Indication**

Orthodontic correction splint and retainer / retention splint / long-term retainer and Mandibular Protrusion Splint / MAD / anti-snoring device and occlusal splint and bruxism splint / miniplast splint

Tips and tricks:

Vacuum-sealed material, does not require predrying. PMMA adjustment can be done without bonding agent, use of ISOFOLAN® required.

Available in:

1.2 / 1.8 / 2.5 / 3.0 / 4.0 mm
x 125 mm Ø

Technical data

Material:	TPU
Hardness:	35 Sh. D
Tensile strength:	39 MPa
Impact strength:	no fracture
Material:	PETG
E-modulus:	2,200 MPa
Hardness:	78 Sh. D
Tensile strength:	53 MPa
Impact strength:	10 KJ/m ²

C € 0044

DURASOFT® blue pd

Indication

Mandibular Protrusion Splint / MAD / anti-snoring device

Tips and tricks:

Vacuum-sealed material, does not require predrying, PMMA adjustment can be done without bonding agent, use of ISOFOLAN® foil required.

Available in:

1.8 / 2.5 mm x 125 mm Ø

Technical data

Material:	TPU
Hardness:	35 Sh. D
Tensile strength:	39 MPa
Impact strength:	no fracture
Material:	PETG
E-modulus:	2,200 MPa
Hardness:	78 Sh. D
Tensile strength:	53 MPa
Impact strength:	10 KJ/m²

DURASOFT® seal

Indication

Discolouration resistant inner layer for DURASOFT® pd splints. Specifically designed for long-term use splints such as Mandibular Protrusion Splint / MAD / anti-snoring device

Tips and tricks:

Can be combined in the lamination process with DURASOFT pd. Use of ISOFOLAN® foil required.

Available in:

0.1 mm x 125 mm Ø
clear, blue, magenta

Technical data

Material:	Thermoplastic polyurethane
E-modulus:	1,900 MPa
Hardness:	82 Sh. D
Tensile strength:	69 MPa
Impact strength:	no fracture

HARDCAST®

Indication

Copings and spacer maintainer for bleaching splints

Tips and tricks:

Copings can be thermoformed with spacer foil.

Available in:

0.4 / 0.6 / 0.8 mm x 125 mm Ø

Technical data

Material:	Polypropylen
E-modulus:	890 MPa
Hardness:	Rockwell 72
Tensile strength:	28 MPa
Impact strength:	600 KJ/m² o. no fracture

IMPRELON® opaque

Indication

Individual impression tray / functional tray

Tips and tricks:

Ready-made handles can be attached using resin, use moistened foam discs as spacer, edges can be smoothed with monomer.

Available in:

2.0 / 3.0 mm x 125 mm Ø
opaque

Technical data

Material:	Polystyrol
E-modulus:	2,500 MPa
Tensile strength:	18 MPa

IMPRELON® clear

Indication

Individual impression tray / functional tray

Tips and tricks:

Ready-made handles can be attached using resin, use moistened foam discs as spacer, edges can be smoothed with monomer.

Available in:

2.0 / 3.0 mm x 125 mm Ø
clear

Technical data

Material:	Polystyrol
E-modulus:	3,200 MPa
Hardness:	Rockwell 150
Tensile strength:	55 MPa
Impact strength:	14 KJ/m²

IMPRELON® S pd / IMPRELON® S+ pd

Indication

Aligner and occlusal splint and retainer / retention splint / Long-term retainer and OSAMU-Retainer®

Tips and tricks:

Can be used in conjunction with the bonding agent OSAMU-Bond, BIOPLAST® 1.5 mm and IMPRELON® S pd 0.75 mm for the production of an OSAMU-Retainer®.

***IMPRELON® S+ pd:** Due to the integrated insulating foil, which is removed after thermoforming, the use of ISOFOLAN® is not required.

Available in:

0.5 / 0.75 / 1.0 / 1.5 / 2.0 / 3.0 mm
x 125 mm Ø

Technical data

Material:	Copolyester
E-modulus:	1,624 MPa
Hardness:	78 Sh. D
Tensile strength:	49 MPa
Impact strength:	no fracture

ISOFOLAN®

Indication

Insulation foil

Tips and tricks:

The ISOFOLAN® foil should be thermoformed, cut and perforated interdentially prior to the fabrication of any kind of hard splints.

Available in:

0.1 mm x 125 mm Ø
natur

Technical data

Material:	HD PE
E-modulus:	600 MPa
Hardness:	45 Sh. D
Tensile strength:	30 MPa
Impact strength:	no fracture

COATING FOIL

Indication

Model protection / protective coating / skin-packaging

Tips and tricks:

Can be used together with the Skin packaging card board for an attractive packaging.

Available in:

0.15 mm x 125 mm Ø

Technical data

Material:	Poly-vinyl-chloride
Hardness:	82 Sh. D
Tensile strength:	46 MPa
Impact strength:	> 550 KJ/m²

SPACER FOIL

Indication

Spacer maintainer foil for copings

Tips and tricks:

Can be used with HARDCAST® or COPYPLAST® for delicate dies.

Available in:

0.1 mm x 125 mm Ø
clear

Technical data

Material:	Poly-vinyl-chloride
Hardness:	78 Sh. D
Tensile strength:	42 MPa
Impact strength:	600 KJ/m ²

UNIQ pd / UNIQ+ pd

Indication

Aligner and retainer / retention splint /
Long-term retainer

Tips and tricks:

Vacuum-sealed foil, no pre-drying
required.

***UNIQ+ pd:** Due to the integrated
insulating foil, which is removed after
thermoforming, the use of ISOFOLAN®
is not required.

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SCHEU-DENTAL GmbH