



# HU-FRIEDY ULTRASONIC TIPS

You trust the premium quality of our scalers and curettes -  
Now experience the same quality with our ultrasonic tips for EMS and Satelec/Acteon generators.

## 2 FOR THE PRICE OF 1 SATISFACTION GUARANTEE

PROMOTION FROM 15<sup>TH</sup> JANUARY TO 30<sup>TH</sup> JUNE 2017

Your satisfaction is at the heart of our business. As such, you can now join thousands of dentists using our Ultrasonic tips for a great price and if you're not delighted with them, you can return them for a full refund.

Example 1: By ordering 1 piece ref. Nr. USS10 and you'll get 2 for the recommended retailprice of 97,25 Euro (48,63 Euro/pc.).

Example 2: By ordering 4 piece ref. Nr. UE100 and you'll get 8 for the recommended retailprice of 488,44 Euro (61,06 Euro/pc.).



Our tips all include a free colour coded torque wrench!



### SERIES S

Compatible with Satelec/Acteon\* generators



**S10 Universal  
USS10**

Standard diameter tip used to remove heavy supragingival deposits.

TOP-SELLER

● POWER SETTING MEDIUM



**100 Universal  
US100**

Probe-sized tip used to remove light deposits in shallow to moderate pockets and for periodontal maintenance.

● POWER SETTING MEDIUM



**P10 Universal  
USP10**

Standard diameter tip used to remove light deposits and plaque. Tip features smooth edges for safe removal of soft and hard deposits.

● POWER SETTING MEDIUM



**3 Perio  
USP3**

Thin-bladed curette-like tip.

● POWER SETTING LOW

### SERIES E

Compatible with all EMS\* generators



**S10 Universal  
UES10**

Standard diameter tip used to remove heavy supragingival deposits.

TOP-SELLER

● POWER SETTING MEDIUM



**100 Universal  
UE100**

Probe-sized tip used to remove light deposits in shallow to moderate pockets and for periodontal maintenance.

● POWER SETTING LOW



**P Universal  
UEP**

Standard diameter tip used to remove moderate to heavy deposits located supragingivally and in shallow pockets.

● POWER SETTING MEDIUM



**3 Perio  
UEP3**

Thin-bladed curette-like tip.

● POWER SETTING LOW

Visit us online at [Hu-Friedy.eu](http://Hu-Friedy.eu)

©2017 Hu-Friedy Mfg. Co., LLC. All rights reserved.

How the best perform




\* Our offered items are compatible with but not original products of device manufacturing companies EMS and/or Satelec and are not produced, licensed or recommended from them. All used brand names and definitions are registered trademarks and brands of respective owners. They are intended only to show compatibility of items with products of different manufacturers.



# HU-FRIEDY ULTRASONIC TIPS

## SERIES S


Compatible with Satelec/Acteon\* generators



Partcode	Description	Reference Satelec	Sugg. Retail-Price	Promo-Price/pc.	Quantity (1+1 package order)
USS10	Piezo Tip S-Series Scaling #S10 universal	S1	97,25 €	48,63 €	
US100	Piezo Tip S-Series Scaling #100 universal	10Z	109,91 €	54,96 €	
USP10	Piezo Tip S-Series Perio #P10 universal	10P	97,25 €	48,63 €	
USP3	Piezo Tip S-Series Perio #3 Perio	H3	115,17 €	57,59 €	

## SERIES E

Compatible with all EMS\* generators



Partcode	Description	Reference EMS	Sugg. Retail-Price	Promo-Price/pc.	Quantity (1+1 package order)
UES10	Piezo Tip E-Series Scaling #S10 universal	A	101,75 €	50,88 €	
UE100	Piezo Tip E-Series Scaling #100 universal	PS	122,11 €	61,06 €	
UEP	Piezo Tip E-Series Scaling #P universal	P	101,75 €	50,88 €	
UEP3	Piezo Tip E-Series Perio #3 Perio		117,02 €	58,51 €	

## GUARDIAN WRENCH/TIP CARRIER

Hu-Friedy Piezo Tips come in the Guardian Wrench Tip Carrier for safe and easy connection to the handpiece. Wrenches also hold the tip during cleaning and sterilisation to prevent loss or damage to the tip.



Scaling Wrench/  
Tip Carrier

### MEDIUM—HIGH

Use primarily on medium to high power (in orange mode or scaling mode if applicable).



Scaling Wrench/  
Tip Carrier

### LOW—MEDIUM

Use primarily on low to medium power (in blue mode or scaling mode if applicable).



Perio Wrench/  
Tip Carrier

### LOW

Use primarily on low power (in green mode or perio mode if applicable).

## HOW IT WORKS

- Complete your name, practice address and contact details along with your preferred Dental Dealer;
- Clearly mark the quantities and items required ( 1 pack = 2 tips);
- On receiving the tips, open just one and test, once satisfied, continue to use all stock;
- In the unlikely event you are not satisfied please send the following document available here ([www.hufriedy.de/hf-returnsform.pdf](http://www.hufriedy.de/hf-returnsform.pdf)), indicate the reason for your dissatisfaction and e-mail to [HFEkundendienst@hu-friedy.com](mailto:HFEkundendienst@hu-friedy.com), return all items to your Dental dealer, please be aware that one tip from each pack MUST be unopened and unused.

Please fill in and fax to Hu-Friedy: Free-Fax 00800 48 37 43 40

Name \_\_\_\_\_

Street \_\_\_\_\_

ZIP/City \_\_\_\_\_

Tel. \_\_\_\_\_

Fax \_\_\_\_\_

E-Mail \_\_\_\_\_

Dental-Depot (Name/City) \_\_\_\_\_

Signature/Stamp \_\_\_\_\_

Please send me regular news about Hu-Friedy products by e-mail

Visit us online at [hu-friedy.eu](http://hu-friedy.eu)

Manufacturer: Hu-Friedy Mfg. Co., LLC | 3232 N. Rockwell Street | Chicago, IL 60618 | USA  
 Customer Care Department: Free Phone 00800 48 37 43 39 | Free Fax 00800 48 37 43 40  
 Website: Hu-Friedy.eu